

# Content of the Tool Set for Installing STable (©) Active Pneumatic Vibration Isolation Tables

Supertech Instruments sends to every new customer a Tool Set for Installing STable (©) Active Pneumatic Vibration Isolation Tables for free only once. It is a free accessory of the first table sold to a reseller or a customer. Further and more Tool Sets can also be ordered. Please see the Price List on our website.

The content of the Tool Set:

4 pieces

Steel rods, made of stainless steel, 250 mm long, 12 mm in diameter

4 pieces

Lifting levers for manual lifting, holding and moving of the top plate, made of steel. They are 200 mm long. There are M16 x 25 mm male screws on one of their ends.

1 piece

8 mm Combined Spanner (Ring and Fork)

1 piece

13 mm Combined Spanner (Ring and Fork)

1 piece

Knife with sharp edge

1 piece

Spirit Level, 200 mm long, very accurate model

1 piece

Hex Socket Wrench Ratchet

1 piece

30 mm Hex Socket Wrench Insert