

# Content of the Tool Set for Installing STable (©) Active Pneumatic Vibration Isolation Tables

Supertech Instruments sends to every new customer a Tool Set for Installing STable (©) Active Pneumatic Vibration Isolation Tables for free only once. It is a free accessory of the first table sold to a reseller or a customer. Further and more Tool Sets can also be ordered. Please see our Price List.

## The content of the Tool Set:

4 pieces

Lifting levers for manual lifting, holding and moving of the top plate, made of steel. They are 150 mm or 200 mm long. There are M16 x 25 mm male threads on one of their ends.

1 piece

8 mm combined wrench (spanner with a ring and a fork end). You should use it to fix the height sensing springs of the valve assemblies after adjustment of the optimal heights.

1 piece

13 mm combined wrench (spanner with a ring and a fork end). You should use it to assemble the arm rest.

1 piece

17 mm combined wrench (spanner with a ring and a fork end). You should use it to adjust the height of the stable feet of the casters.

1 piece

PH2 screwdriver. You should use it to remove, then reassemble the stainless steel covers of the front valve assemblies.

1 piece

Spirit Level, 200 mm long, very accurate model. You should use it to measure the horizontal position of the lower horizontal bars of the leg system. You can adjust the horizontality of the leg system by the stable feet of the casters.

For the details of the adjustments please read through carefully both the User Manual and the Installation Manual. The installation and the adjustment of the table are not difficult. But there are a few points in the process, which you must follow carefully if you do not want to make **serious damage to the table**.