

## **Elastomer-based Damping Plates**

Affordable Vibration Isolation



### **General Description:**

Damping Plates are elastomer-based vibration isolator plates designed for applications that do not require high level of vibration attenuation. They can be put onto any given table. Our Damping Plates are available in *two fixed sizes*:

**60 x 40 cm**

Height: 90 mm

Mass load: 300 N

**75 x 50 cm**

Height: 90 mm

Mass load: 500 N

### **Specifications:**

- painted carbon steel top plate, magnetic
- vibration isolation is passive
- vibrations are absorbed by elastomer (silicone rubber) cones
- plate can be leveled by adjusting the legs
- plate has a build-in spirit level so leveling can be corrected even if the load is not concentric



### **Warranty:**

Supertech Ltd gives you 3 years of full warranty for its mechanical products against defects in materials or workmanship as long as the equipment has been subjected to normal and proper use. During this period, faulty products will be repaired or replaced free of charge provided they are returned to our workshop. (Postage prepaid by the Customer.) Supertech Ltd will undertake the servicing and calibration after the expiration of the warranty period for a nominal fee.