



Triangle:
Sign of PIN(1)



- 2 PIN(1): Out GND
- PIN(2): Output 1, or 9
- PIN(3): Out GND
- 4 PIN(4): Output 2, or 10
- PIN(5): Out GND
- 6 PIN(6): Output 3, or 11
- PIN(7): Out GND
- 8 PIN(8): Output 4, or 12
- PIN(9): Out GND
- 10 PIN(10): Output 5, or 13
- PIN(11): Out GND
- 12 PIN(12): Output 6, or 14
- PIN(13): Out GND
- 14 PIN(14): Output 7, or 15
- PIN(15): Out GND
- 16 PIN(16): Output 8, or 16
- PIN(17): Out GND
- PIN(18): Out GND
- 18 PIN(19): Out GND
- PIN(20): Out GND

The appropriate 20-pin, 2.54 mm female plug compatible with this output connector is IDS-020-S100-01/P manufactured by E-tec Corp.

Output Connector of each 8-channel sections of MultiAmp, V.2